

Adult 7" Length Transtibial APOPPS® Measurement Form

DATE: _____

PATIENT: _____ PROSTHETIST: _____ PO NUMBER: _____

MEASUREMENTS (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM: Left** _____ **Right** _____

LENGTH: _____ In. (Cm) MPT TO DISTAL END (shorter than 4" or longer than 8 1/4" without bandages requires custom made sockets)

CIRCUMFERENCES: _____ In. (Cm) MPT _____ In. (Cm) 2" UP FROM DISTAL END (Choose size based upon distal circ.)

Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered. CUSTOM _____

Use the following sizing information for selecting APOPPS™ sockets for your patient. Each of these sockets is a one piece, polyethylene, clam shell design with anatomically trimmed overlapping sides, ideal for: accommodating bulbous, straight or conical shapes, protecting and gently reducing the residuum throughout the rehabilitation period and for weight bearing when used in conjunction with the UFOS™.

- The FLO-TECH-TOR™ is used as a protective socket and for early supervised weight bearing prosthetic management of transtibial amputations.
- The VCSPS™ in conjunction with the UFOS™ is used for preparatory prosthetic management of transtibial amputations.
- The UFOS™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR™ or VCSPS™ socketed. This system allows your patient to wear the FLO-TECH-TOR™ or VCSPS™ sockets alone when the pylon and foot are not needed.

The APOPPS™ System is designed to keep a healthy amputation healthy throughout the rehabilitation process.

ADULT SIZES: *FLO-TECH-TOR™ VCSPS™ UFOS™ REHAB PREP *APOPPS®

7" (4 1/4" - 6 1/4") **LENGTH**

X-Small: MPT 14" to 15" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 13" (33cm) to 14" (36cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 11" (28cm) to 15" (38cm)

Small: MPT 15" to 16" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 14" (36cm) to 15" (38cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 12" (30cm) to 16" (41cm)

Medium: MPT 16" to 17" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 15" (38cm) to 16" (41cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 13" (33cm) to 17" (43cm)

Large: MPT 17" to 18" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 16" (41cm) to 17" (43cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 14" (36cm) to 18" (46cm)

X-Large: MPT 18" to 19" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 17" (43cm) to 18" (46cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 15" (38cm) to 20" (51cm)

XX-Large: MPT 19" to 20" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 18" (46cm) to 19" (48cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 16" (41cm) to 22" (53cm)

XXX-Large: MPT 21" to 22" 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____ 7-L _____

Distal: 20" (51cm) to 21" (53.5cm) 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____ 7-R _____

Range: 19" (48.5cm) to 23" (58.5cm)

*DOES NOT INCLUDE RETICULATED PAD

(9" length available with separate form, ALL OTHER DIMENSIONS & LENGTHS ARE CUSTOM ORDERS)

COMPLETE FOOT & PYLON SYSTEMS

Complete Titanium Pylon System _____

Complete Stainless Steel Pylon System _____

(Required with all foot orders)

Includes: Pyramid Plate, S.A.C.H. Foot Adapter, 30mm Pylon, 30mm Tube Clamp **AND S.A.C.H. Foot:** L R **Size:** 24 26 28

COMFORT OPTIONS & COMPONENT ITEMS:

QTY	DESCRIPTION
_____	Modification to Standard Socket
_____	Waist Belt/Suspension Strap
_____	Fork/Suspension Strap
_____	Per Socket Strap - ALL SIZES
_____	Velcro® Compatible Neoprene® socket band _____
_____	Distal End Pads - Per pkg. SIZE: _____
_____	Distal End Pad, Reticulated SIZE: _____
_____	PowerPoint Marketing Presentation
_____	Marketing Packets (10 pack)
_____	CONSUMER GUIDES (10 pack)
_____	NEW! MARKETING KIT
_____	Waist Belt/Suspension Strap
_____	Catalog (one per facility)

QTY	DESCRIPTION
_____	Each Amputation Socks SIZE: _____
_____	4 Bolt Adapter Plate, 4 short, 4 long screws
_____	Post-op Foot - L _____, R _____, SIZE: _____
_____	Titanium SACH Foot Adapter
_____	Titanium Pyramid Skt. Attachment Plate
_____	Titanium 30 mm Pylon Tube _____
_____	Titanium 30 mm Tube Clamp
_____	SS SACH Foot Adapter
_____	SS Pyramid Socket attachment Plate
_____	SS 30 mm Pylon Tube
_____	SS 30 mm Tube Clamp
_____	Pour & Modify Cast
_____	CONSUMER GUIDE (each)



FLO-TECH® O & P Systems, Inc.
 PO Box 462 - 7325 Halseyville Rd
 Trumansburg, New York 14886

www.1800FLO-TECH.com
info@1800FLO-TECH.com

Phone: (607) 387-3070

Phone: 1-800-FLO-TECH (356-8324)

Fax: (607) 387-3176

Adult 9" Length Transtibial APOPPS® Measurement Form

DATE: _____

PATIENT: _____ PROSTHETIST: _____ PO NUMBER: _____

MEASUREMENTS (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM: Left** _____ **Right** _____

LENGTH: _____ In. (Cm) MPT TO DISTAL END (shorter than 4" or longer than 8 1/4" without bandages requires custom made sockets)

CIRCUMFERENCES: _____ In. (Cm) MPT _____ In. (Cm) 2" UP FROM DISTAL END (Choose size based upon distal circ.)

Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered. CUSTOM _____

Use the following sizing information for selecting APOPPS™ sockets for your patient. Each of these sockets is a one piece, polyethylene, clam shell design with anatomically trimmed overlapping sides, ideal for: accommodating bulbous, straight or conical shapes, protecting and gently reducing the residuum throughout the rehabilitation period and for weight bearing when used in conjunction with the UFOS™.

- The FLO-TECH-TOR™ is used as a protective socket and for early supervised weight bearing prosthetic management of transtibial amputations.
- The VCSPS™ in conjunction with the UFOS™ is used for preparatory prosthetic management of transtibial amputations.
- The UFOS™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR™ or VCSPS™ socketed. This system allows your patient to wear the FLO-TECH-TOR™ or VCSPS™ sockets alone when the pylon and foot are not needed.

The APOPPS™ System is designed to keep a healthy amputation healthy throughout the rehabilitation process.

ADULT SIZES: *FLO-TECH-TOR™ VCSPS™ UFOS™ REHAB PREP *APOPPS®

9" (6 1/2" - 8 1/4") Length

X-Small: MPT 14" to 15" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 13" (33cm) to 14" (36cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 11" (28cm) to 15" (38cm)

Small: MPT 15" to 16" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 14" (36cm) to 15" (38cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 12" (30cm) to 16" (41cm)

Medium: MPT 16" to 17" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 15" (38cm) to 16" (41cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 13" (33cm) to 17" (43cm)

Large: MPT 17" to 18" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 16" (41cm) to 17" (43cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 14" (36cm) to 18" (46cm)

X-Large: MPT 18" to 19" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 17" (43cm) to 18" (46cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 15" (38cm) to 20" (51cm)

XX-Large: MPT 19" to 20" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 18" (46cm) to 19" (48cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 16" (41cm) to 22" (53cm)

XXX-Large: MPT 21" to 22" 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____ 9-L _____

Distal: 20" (51cm) to 21" (53.5cm) 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____ 9-R _____

Range: 19" (48.5cm) to 23" (58.5cm)

*DOES NOT INCLUDE RETICULATED PAD (7" length available with separate form, ALL OTHER DIMENSIONS & LENGTHS ARE CUSTOM ORDERS)

COMPLETE FOOT & PYLON SYSTEMS Complete Titanium Pylon System _____ Complete Stainless Steel Pylon System _____ (Required with all foot orders)

Includes: Pyramid Plate, S.A.C.H. Foot Adapter, 30mm Pylon, 30mm Tube Clamp **AND S.A.C.H. Foot:** L R **Size:** 24 26 28

COMFORT OPTIONS & COMPONENT ITEMS:

QTY DESCRIPTION	QTY DESCRIPTION
_____ Modification to Standard Socket	_____ Each Amputation Socks SIZE: _____
_____ Waist Belt/Suspension Strap	_____ 4 Bolt Adapter Plate, 4 short, 4 long screws
_____ Fork/Suspension Strap	_____ Post-op Foot - L ____, R ____, SIZE: _____
_____ Per Socket Strap - ALL SIZES	_____ Titanium SACH Foot Adapter
_____ Velcro® Compatible Neoprene® socket band _____	_____ Titanium Pyramid Skt. Attachment Plate
_____ Distal End Pads - Per pkg. SIZE: _____	_____ Titanium 30 mm Pylon Tube _____
_____ Distal End Pad, Reticulated SIZE: _____	_____ Titanium 30 mm Tube Clamp
_____ PowerPoint Marketing Presentation	_____ SS SACH Foot Adapter
_____ Marketing Packets (10 pack)	_____ SS Pyramid Socket attachment Plate
_____ CONSUMER GUIDES (10 pack)	_____ SS 30 mm Pylon Tube
_____ NEW! MARKETING KIT	_____ SS 30 mm Tube Clamp
_____ Waist Belt/Suspension Strap	_____ Pour & Modify Cast
_____ Catalog (one per facility)	_____ CONSUMER GUIDE (each)



FLO-TECH® O & P Systems, Inc.
PO Box 462 - 7325 Halseyville Rd
Trumansburg, New York 14886

www.1800FLO-TECH.com
info@1800FLO-TECH.com

Phone: (607) 387-3070

Phone: 1-800-FLO-TECH (356-8324)

Fax: (607) 387-3176

Transfemoral APOPPS-TF® Measurement Form

DATE: _____

PATIENT: _____ PROSTHETIST: _____ PO NUMBER: _____

Each of these sockets is a one piece, polyethylene, clam shell design with overlapping sides, ideal for accommodating bulbous, straight or conical shapes and for gently reducing the residuum throughout the rehabilitation period.

- The FLO-TECH-TOR-TF™ is used as a protective socket and for early supervised weight bearing prosthetic management of transfemoral amputations when used in conjunction with the UFOS™.
- The UFOS-TF™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR-TF™ you select for your patient. This system allows your patient to wear the FLO-TECH-TOR-TF™ alone when the pylon and foot are not needed.

MEASUREMENTS (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.) _____ **CUSTOM**

LENGTH: (ISCHIUM TO DISTAL END OVER BANDAGES) _____ In. (Cm), **RESIDUUM:** _____ **Left** _____ **Right**

***CIRCUMFERENCE** (ISCHIAL LEVEL) _____ In. (Cm), (3" UP FROM DISTAL END) _____ In. (Cm)

***SOCKETS ARE SIZED TO FIT IDEAL CIRCUMFERENCE AT THE ABOVE ISCHIAL LOCATION. (Top measurement at each size indicates ideal Ischial level circumference; bottom measurement indicates the minimum to maximum circumference range at the Ischial level.)**

*SOCKETS ARE SIZED TO FIT THIS RANGE AT THE ISCHIAL LEVEL MAXIMUM LENGTH 11 1/4" – THEN CUSTOM

**SOCKETS ARE SIZED TO FIT THIS RANGE AT THE DISTAL END

***DOES NOT INCLUDE RETICULATED PAD

SOCKET AND/OR PACKAGE:

12" (9 3/4" - 11 1/4") Length

***FLO-TECH-TOR-TF™

UFOS-TF™

***APOPPS-TF®

SMALL: *16" (41cm) to 22" (56cm)

L _____

L _____

L _____

** Range 13" (33cm) to 17" (31cm)

R _____

R _____

R _____

MEDIUM: *19" (48cm) to 25" (64cm)

L _____

L _____

L _____

** Range 16" (40cm) to 20" (51cm)

R _____

R _____

R _____

LARGE: *22" (56cm) to 28" (71cm)

L _____

L _____

L _____

** Range 19" (48cm) to 23" (58cm)

R _____

R _____

R _____

X-LARGE: *26" (66cm) to 32" (81cm)

L _____

L _____

L _____

** Range 24" (61cm) to 28" (71cm)

R _____

R _____

R _____

(All Other Dimensions and Lengths Are Custom Orders)

L-CODES (ADD ALL OTHERS YOU PROVIDE) L-5460 OR L-5420 OR L-5585

SYSTEM OPTIONS:

Complete Ti/Alum. Knee/Pylon System _____

Complete SS/Alum. Knee/Pylon System _____

Includes: Locking Knee, Rotatable 4-hole Connector, S.A.C.H. Foot Adapter, 40mm Pylon **AND**

S.A.C.H. Foot: L R **Size:** 24 26 28 (Required with all foot orders)

COMFORT OPTIONS & COMPONENT ITEMS:

QTY DESCRIPTION

- ____ Modification to Standard Socket
- ____ Neoplush Pelvic Suspension band (2-Pc.)
- ____ Neoplush Pelvic Band (Amputated Side)
- ____ Neoplush Pelvic Band (Contra-lateral Side)
- ____ Velcro® Compatible Neoprene® Bands - ALL TYPES
- ____ Per Socket Strap - ALL SIZES
- ____ Distal End Pad Kit (1" & 2") SIZE: _____
- ____ Distal End Pad, Reticulated SIZE: _____
- ____ **NEW!** Marketing Kit
- ____ PowerPoint Marketing Presentation
- ____ Marketing Packets (10 pack)
- ____ Selsian Bandage [Belt – no strap]
- ____ Selsian Bandage w/Auxiliary Strap
- ____ Catalog (one per facility)

QTY DESCRIPTION

- ____ Each Amputation Socks SIZE: _____
- ____ 4 Bolt Adapter Plate, 4 short, 4 long screws
- ____ Post-op Foot - L R SIZE: _____
- ____ Titanium S.A.C.H. Foot Adapter
- ____ Titanium Rotatable 4-hole Connector
- ____ Titanium 30mm Pylon Tube
- ____ Aluminum Manual Locking Knee (<300#)
- ____ SS S.A.C.H. Foot Adapter
- ____ SS Rotatable 4-hole Connector
- ____ SS 30 mm Pylon Tube
- ____ Pour & Modify Cast
- ____ CONSUMER GUIDE (each)
- ____ CONSUMER GUIDES (10 pack)



FLO-TECH® O & P Systems, Inc.
 PO Box 462 - 7325 Halseyville Rd
 Trumansburg, New York 14886

www.1800FLO-TECH.com
info@1800FLO-TECH.com

Phone: (607) 387-3070

Phone: 1-800-FLO-TECH (356-8324)

Fax: (607) 387-3176

Measurement Form

VISA/MASTERCARD ACCT. # _____ EXP. DATE: _____

NAME AS IT APPEARS ON CARD: _____ INVOICE#: _____

PATIENT: _____ PROSTHETIST: _____

COMPANY: _____ NEED BY: _____

PHONE: _____ FAX: _____

Shipping Address: _____ Billing Address: _____

TRANSTIBIAL

10" Length (4½" – 10") **ONE SIZE** **LEFT _____ RIGHT _____**

Distal Range 11" (28cm) to 23" (58½cm)

1. A modified FLO-TECH-TOR™ made from 3/16" firm Pelite
 - a. Hook Velcro®
2. A standard FLO-TECH® Neoprene Thigh Band
3. An auxiliary Neoprene Band applied at the knee
4. A standard 9" FLO-TECH® Neoprene Socket Band
5. A Neoprene Distal Socket Retainer (cap)
6. Standard 1" distal end pad
7. Standard 2" distal end pad

TRANSFEMORAL

13" Length (7½" - 13") **ONE SIZE** **UNIVERSAL SIDE**
Distal Range 13" (33cm) to 23" (58½cm) **QUANTITY** _____

1. A FLO-TECH-TOR-TF™ made from ¾" firm Pelite
 - a. Hook Velcro®
2. A standard FLO-TECH® Neoprene 2-piece Pelvic Belt
3. A 9" Lateral Socket Extension Band (used to extend the lateral Trochanter area of the socket to interface with the Trochanter section of the Pelvic Belt)
 - a. A 2" x 5" double sided hook Velcro® pressed between a standard Neoprene Trochanter section of the pelvic belt and the Lateral Socket Extension Band
 - b. The Lateral Socket Extension Band may be used horizontally or vertically depending on the size of the patient
4. It is attached to the SOFTIE-TF™ via 2 patches of Velcro® at the Proximal Lateral area
5. A standard FLO-TECH-TOR-TF™ Neoprene Socket Band
6. A Neoprene Distal Socket Retainer
7. (2) Standard 1" distal end pads
8. (2) Standard 2" distal end pads

Available Options:

QTY DESCRIPTION	QTY DESCRIPTION
____ Transtibial waistbelt and forkstrap	____ Amputation Sock SIZE: _____
____ Suspension strap	____ CATALOG (One Per Facility)
____ Neoplush pelvic band	____ AMPUTEE GUIDE - (Each)
____ Neoplush contra-lateral band	____ AMPUTEE GUIDES - (10 Pack)
____ Neoplush socket band	____ Marketing Packets - (10 Pack)
____ Distal end pads/Per Pkg.	____ NEW! MARKETING KIT
____ Distal end pad/Reticulated SIZE: _____	____ PowerPoint Marketing Presentation



FLO-TECH® O & P Systems, Inc.
 PO Box 462 - 7325 Halseyville Rd
 Trumansburg, New York 14886

www.1800FLO-TECH.com
info@1800FLO-TECH.com
Phone: (607) 387-3070
Phone: 1-800-FLO-TECH(356-8324)
Fax: (607) 387-3176

SYMES Measurement Form

DATE: _____ P.O. Number: _____

VISA/MASTERCARD ACCT. # _____ EXP. DATE: _____

NAME AS IT APPEARS ON CARD: _____ INVOICE#: _____

PATIENT: _____ PROSTHETIST: _____

COMPANY: _____ NEED BY: _____

PHONE: _____ FAX: _____

Shipping Address: _____ Billing Address: _____

CUSTOM SYMES ORDER: (PROVIDED BY PRACTITIONER)

LEFT: _____ RIGHT: _____ From Measurements: _____ From Tracing: _____ (FAX Two 8" x 11" pieces of paper)

CIRCLE SOCKETS FOR THIS ORDER:

LENGTHS:

From Mid Patella Tendon To Distal End _____ IN _____ CM

From Smallest Part of Calf to Distal End _____ IN _____ CM

From Smallest Part Ankle to Distal End _____ IN _____ CM

AP (ANTERIOR TO POSTERIOR):

At Smallest Part of Ankle _____ IN _____ CM

At Ball of Symes _____ IN _____ CM

ML (MEDIAL TO LATERAL SIDE):

At Tibial Tubercle _____ IN _____ CM

At Largest Part of Calf _____ IN _____ CM

At Smallest Part of Ankle _____ IN _____ CM

At Ball of Symes _____ IN _____ CM

CIRCUMFERENCES:

At Tibia Tubercle _____ IN _____ CM

At Widest Part of Calf _____ IN _____ CM

At Just below Calf Muscle _____ IN _____ CM

At Mid Point From Distal Calf to Smallest Part of Ankle _____ IN _____ CM

At Smallest Part of Ankle _____ IN _____ CM

At Ball of Symes amputation _____ IN _____ CM

COMFORT OPTIONS & COMPONENT ITEM PRICES:

QTY	DESCRIPTION
_____	Modification to Standard Socket
_____	Distal End Pads - Per pkg. SIZE: _____
_____	Waist Belt/Suspension Strap
_____	Pour & Modify Cast
_____	Fork/Suspension Strap
_____	Per Socket Strap - ALL SIZES
_____	PowerPoint Marketing Presentation
_____	Velcro® Compatible Neoprene® socket band _____
_____	NEW! MARKETING KIT

QTY	DESCRIPTION
_____	Each Amputation Sock SIZE: _____
_____	SS Rotatable Connector
_____	SS SACH Foot Adapter
_____	Foot - L _____, R _____, SIZE _____
_____	Titanium Rotatable 4-hole Connector
_____	Titanium SACH Foot Adapter
_____	Titanium SACH Foot Adapter
_____	Marketing Packets (10 pack)
_____	CONSUMER GUIDES (10 pack)



FLO-TECH® O & P Systems, Inc.
 PO Box 462 - 7325 Halseyville Rd
 Trumansburg, New York 14886
 www.1800FLO-TECH.com

All materials © FLO-TECH® May 11, 1995, Modified July 2008

Phone: 1-800-FLO-TECH (356-8324)
Phone: (607) 387-3070
Fax: (607) 387-3176
E-mail: info@1800FLO-TECH.com