



# Adult 9" Length Transtibial APOPPS® Measurement Form

DATE: \_\_\_\_\_

PATIENT: \_\_\_\_\_ PROSTHETIST: \_\_\_\_\_ PO NUMBER: \_\_\_\_\_

**MEASUREMENTS** (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM: Left** \_\_\_\_\_ **Right** \_\_\_\_\_

**LENGTH:** \_\_\_\_\_ In. (Cm) MPT TO DISTAL END (shorter than 4" or longer than 8 1/4" without bandages requires custom made sockets)

**CIRCUMFERENCES:** \_\_\_\_\_ In. (Cm) MPT \_\_\_\_\_ In. (Cm) 2" UP FROM DISTAL END (Choose size based upon distal circ.)

Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered. CUSTOM \_\_\_\_\_

Use the following sizing information for selecting APOPPS™ sockets for your patient. Each of these sockets is a one piece, polyethylene, clam shell design with anatomically trimmed overlapping sides, ideal for: accommodating bulbous, straight or conical shapes, protecting and gently reducing the residuum throughout the rehabilitation period and for weight bearing when used in conjunction with the UFOS™.

- The FLO-TECH-TOR™ is used as a protective socket and for early supervised weight bearing prosthetic management of transtibial amputations.
- The VCSPS™ in conjunction with the UFOS™ is used for preparatory prosthetic management of transtibial amputations.
- The UFOS™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR™ or VCSPS™ socketed. This system allows your patient to wear the FLO-TECH-TOR™ or VCSPS™ sockets alone when the pylon and foot are not needed.

*The APOPPS™ System is designed to keep a healthy amputation healthy throughout the rehabilitation process.*

**ADULT SIZES:** \_\_\_\_\_ \*FLO-TECH-TOR™ VCSPS™ UFOS™ REHAB PREP \*APOPPS®

**9" (6 1/2" - 8 1/4") Length**

**X-Small:** MPT 14" to 15" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 13" (33cm) to 14" (36cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 11" (28cm) to 15" (38cm)

**Small:** MPT 15" to 16" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 14" (36cm) to 15" (38cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 12" (30cm) to 16" (41cm)

**Medium:** MPT 16" to 17" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 15" (38cm) to 16" (41cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 13" (33cm) to 17" (43cm)

**Large:** MPT 17" to 18" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 16" (41cm) to 17" (43cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 14" (36cm) to 18" (46cm)

**X-Large:** MPT 18" to 19" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 17" (43cm) to 18" (46cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 15" (38cm) to 20" (51cm)

**XX-Large:** MPT 19" to 20" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 18" (46cm) to 19" (48cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 16" (41cm) to 22" (53cm)

**XXX-Large:** MPT 21" to 22" 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_ 9-L \_\_\_\_\_

Distal: 20" (51cm) to 21" (53.5cm) 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_ 9-R \_\_\_\_\_

Range: 19" (48.5cm) to 23" (58.5cm)

\*DOES NOT INCLUDE RETICULATED PAD

(7" length available with separate form, ALL OTHER DIMENSIONS & LENGTHS ARE CUSTOM ORDERS)

## COMPLETE FOOT & PYLON SYSTEMS

Complete Titanium Pylon System \_\_\_\_\_

Complete Stainless Steel Pylon System \_\_\_\_\_

(Required with all foot orders)

**Includes:** Pyramid Plate, S.A.C.H. Foot Adapter, 30mm Pylon, 30mm Tube Clamp **AND S.A.C.H. Foot:**  L  R **Size:**  24  26  28

## COMFORT OPTIONS & COMPONENT ITEMS:

- | QTY   | DESCRIPTION                                    |
|-------|--|
| _____ | Modification to Standard Socket                |
| _____ | Waist Belt/Suspension Strap                    |
| _____ | Fork/Suspension Strap                          |
| _____ | Per Socket Strap - ALL SIZES                   |
| _____ | Velcro® Compatible Neoprene® socket band _____ |
| _____ | Distal End Pads - Per pkg. SIZE: _____         |
| _____ | Distal End Pad, Reticulated SIZE: _____        |
| _____ | PowerPoint Marketing Presentation              |
| _____ | Marketing Packets (10 pack)                    |
| _____ | CONSUMER GUIDES (10 pack)                      |
| _____ | NEW! MARKETING KIT                             |
| _____ | Waist Belt/Suspension Strap                    |
| _____ | Catalog (one per facility)                     |

- | QTY   | DESCRIPTION                                  |
|-------|--|
| _____ | Each Amputation Socks SIZE: _____            |
| _____ | 4 Bolt Adapter Plate, 4 short, 4 long screws |
| _____ | Post-op Foot - L ____, R ____, SIZE: _____   |
| _____ | Titanium SACH Foot Adapter                   |
| _____ | Titanium Pyramid Skt. Attachment Plate       |
| _____ | Titanium 30 mm Pylon Tube _____              |
| _____ | Titanium 30 mm Tube Clamp                    |
| _____ | SS SACH Foot Adapter                         |
| _____ | SS Pyramid Socket attachment Plate           |
| _____ | SS 30 mm Pylon Tube                          |
| _____ | SS 30 mm Tube Clamp                          |
| _____ | Pour & Modify Cast                           |
| _____ | CONSUMER GUIDE (each)                        |



**FLO-TECH® O & P Systems, Inc.**  
 PO Box 462 - 7325 Halseyville Rd  
 Trumansburg, New York 14886

[www.1800FLO-TECH.com](http://www.1800FLO-TECH.com)  
[info@1800FLO-TECH.com](mailto:info@1800FLO-TECH.com)

Phone: (607) 387-3070

Phone: 1-800-FLO-TECH (356-8324)

Fax: (607) 387-3176

# Transfemoral APOPPS-TF® Measurement Form

DATE: \_\_\_\_\_

PATIENT: \_\_\_\_\_ PROSTHETIST: \_\_\_\_\_ PO NUMBER: \_\_\_\_\_

Each of these sockets is a one piece, polyethylene, clam shell design with overlapping sides, ideal for accommodating bulbous, straight or conical shapes and for gently reducing the residuum throughout the rehabilitation period.

- The FLO-TECH-TOR-TF™ is used as a protective socket and for early supervised weight bearing prosthetic management of transfemoral amputations when used in conjunction with the UFOS™.
- The UFOS-TF™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR-TF™ you select for your patient. This system allows your patient to wear the FLO-TECH-TOR-TF™ alone when the pylon and foot are not needed.

**MEASUREMENTS** (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.) \_\_\_\_\_ **CUSTOM**

**LENGTH:** (ISCHIUM TO DISTAL END OVER BANDAGES) \_\_\_\_\_ In. (Cm), **RESIDUUM:** \_\_\_\_\_ **Left** \_\_\_\_\_ **Right**

\***CIRCUMFERENCE** (ISCHIAL LEVEL) \_\_\_\_\_ In. (Cm), (3" UP FROM DISTAL END) \_\_\_\_\_ In. (Cm)

**\*SOCKETS ARE SIZED TO FIT IDEAL CIRCUMFERENCE AT THE ABOVE ISCHIAL LOCATION. (Top measurement at each size indicates ideal Ischial level circumference; bottom measurement indicates the minimum to maximum circumference range at the Ischial level.)**

\*SOCKETS ARE SIZED TO FIT THIS RANGE AT THE ISCHIAL LEVEL MAXIMUM LENGTH 11 1/4" – THEN CUSTOM

\*\*SOCKETS ARE SIZED TO FIT THIS RANGE AT THE DISTAL END

\*\*\*DOES NOT INCLUDE RETICULATED PAD

**SOCKET AND/OR PACKAGE:**

**12" (9 3/4" - 11 1/4") Length**

\*\*\*FLO-TECH-TOR-TF™

UFOS-TF™

\*\*\*APOPPS-TF®

**SMALL:** \*16" (41cm) to 22" (56cm)

L \_\_\_\_\_

L \_\_\_\_\_

L \_\_\_\_\_

\*\* Range 13" (33cm) to 17" (31cm)

R \_\_\_\_\_

R \_\_\_\_\_

R \_\_\_\_\_

**MEDIUM:** \*19" (48cm) to 25" (64cm)

L \_\_\_\_\_

L \_\_\_\_\_

L \_\_\_\_\_

\*\* Range 16" (40cm) to 20" (51cm)

R \_\_\_\_\_

R \_\_\_\_\_

R \_\_\_\_\_

**LARGE:** \*22" (56cm) to 28" (71cm)

L \_\_\_\_\_

L \_\_\_\_\_

L \_\_\_\_\_

\*\* Range 19" (48cm) to 23" (58cm)

R \_\_\_\_\_

R \_\_\_\_\_

R \_\_\_\_\_

**X-LARGE:** \*26" (66cm) to 32" (81cm)

L \_\_\_\_\_

L \_\_\_\_\_

L \_\_\_\_\_

\*\* Range 24" (61cm) to 28" (71cm)

R \_\_\_\_\_

R \_\_\_\_\_

R \_\_\_\_\_

(All Other Dimensions and Lengths Are Custom Orders)

**L-CODES** (ADD ALL OTHERS YOU PROVIDE) L-5460 OR L-5420 OR L-5585

**SYSTEM OPTIONS:**

Complete Ti/Alum. Knee/Pylon System \_\_\_\_\_

Complete SS/Alum. Knee/Pylon System \_\_\_\_\_

**Includes:** Locking Knee, Rotatable 4-hole Connector, S.A.C.H. Foot Adapter, 40mm Pylon **AND**

**S.A.C.H. Foot:**  L  R **Size:**  24  26  28 (Required with all foot orders)

**COMFORT OPTIONS & COMPONENT ITEMS:**

**QTY DESCRIPTION**

- \_\_\_\_ Modification to Standard Socket
- \_\_\_\_ Neoplush Pelvic Suspension band (2-Pc.)
- \_\_\_\_ Neoplush Pelvic Band (Amputated Side)
- \_\_\_\_ Neoplush Pelvic Band (Contra-lateral Side)
- \_\_\_\_ Velcro® Compatible Neoprene® Bands - ALL TYPES
- \_\_\_\_ Per Socket Strap - ALL SIZES
- \_\_\_\_ Distal End Pad Kit (1" & 2") SIZE: \_\_\_\_\_
- \_\_\_\_ Distal End Pad, Reticulated SIZE: \_\_\_\_\_
- \_\_\_\_ **NEW!** Marketing Kit
- \_\_\_\_ PowerPoint Marketing Presentation
- \_\_\_\_ Marketing Packets (10 pack)
- \_\_\_\_ Selsian Bandage [Belt – no strap]
- \_\_\_\_ Selsian Bandage w/Auxiliary Strap
- \_\_\_\_ Catalog (one per facility)

**QTY DESCRIPTION**

- \_\_\_\_ Each Amputation Socks SIZE: \_\_\_\_\_
- \_\_\_\_ 4 Bolt Adapter Plate, 4 short, 4 long screws
- \_\_\_\_ Post-op Foot - L  R  SIZE: \_\_\_\_\_
- \_\_\_\_ Titanium S.A.C.H. Foot Adapter
- \_\_\_\_ Titanium Rotatable 4-hole Connector
- \_\_\_\_ Titanium 30mm Pylon Tube
- \_\_\_\_ Aluminum Manual Locking Knee (<300#)
- \_\_\_\_ SS S.A.C.H. Foot Adapter
- \_\_\_\_ SS Rotatable 4-hole Connector
- \_\_\_\_ SS 30 mm Pylon Tube
- \_\_\_\_ Pour & Modify Cast
- \_\_\_\_ CONSUMER GUIDE (each)
- \_\_\_\_ CONSUMER GUIDES (10 pack)



**FLO-TECH® O & P Systems, Inc.**  
 PO Box 462 - 7325 Halseyville Rd  
 Trumansburg, New York 14886

**www.1800FLO-TECH.com**  
**info@1800FLO-TECH.com**

**Phone: (607) 387-3070**

**Phone: 1-800-FLO-TECH (356-8324)**

**Fax: (607) 387-3176**

**SOFTIE™ & SOFTIE-TF™**

DATE: \_\_\_\_\_ P.O. Number: \_\_\_\_\_

**Measurement Form**

VISA/MASTERCARD ACCT. # \_\_\_\_\_ EXP. DATE: \_\_\_\_\_

NAME AS IT APPEARS ON CARD: \_\_\_\_\_ INVOICE#: \_\_\_\_\_

PATIENT: \_\_\_\_\_ PROSTHETIST: \_\_\_\_\_

COMPANY: \_\_\_\_\_ NEED BY: \_\_\_\_\_

PHONE: \_\_\_\_\_ FAX: \_\_\_\_\_

Shipping Address: \_\_\_\_\_ Billing Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**TRANSTIBIAL**

**10" Length (4½" – 10")** **ONE SIZE** **LEFT \_\_\_\_\_ RIGHT \_\_\_\_\_**

**Distal Range 11" (28cm) to 23" (58½cm)**

1. A modified FLO-TECH-TOR™ made from 3/16" firm Pelite
  - a. Hook Velcro®
2. A standard FLO-TECH® Neoprene Thigh Band
3. An auxiliary Neoprene Band applied at the knee
4. A standard 9" FLO-TECH® Neoprene Socket Band
5. A Neoprene Distal Socket Retainer (cap)
6. Standard 1" distal end pad
7. Standard 2" distal end pad

**TRANSFEMORAL**

**13" Length (7½" - 13")** **ONE SIZE** **UNIVERSAL SIDE**

**Distal Range 13" (33cm) to 23" (58½cm)**

1. A FLO-TECH-TOR-TF™ made from ¾" firm Pelite
  - a. Hook Velcro®
2. A standard FLO-TECH® Neoprene 2-piece Pelvic Belt
3. A 9" Lateral Socket Extension Band (used to extend the lateral Trochanter area of the socket to interface with the Trochanter section of the Pelvic Belt)
  - a. A 2" x 5" double sided hook Velcro® pressed between a standard Neoprene Trochanter section of the pelvic belt and the Lateral Socket Extension Band
  - b. The Lateral Socket Extension Band may be used horizontally or vertically depending on the size of the patient
4. It is attached to the SOFTIE-TF™ via 2 patches of Velcro® at the Proximal Lateral area
5. A standard FLO-TECH-TOR-TF™ Neoprene Socket Band
6. A Neoprene Distal Socket Retainer
7. (2) Standard 1" distal end pads
8. (2) Standard 2" distal end pads

**Available Options:**

| QTY DESCRIPTION                             | QTY DESCRIPTION                        |
|---|--|
| ____ Transtibial waistbelt and forkstrap    | ____ Amputation Sock SIZE: _____       |
| ____ Suspension strap                       | ____ CATALOG (One Per Facility)        |
| ____ Neoplush pelvic band                   | ____ AMPUTEE GUIDE - (Each)            |
| ____ Neoplush contra-lateral band           | ____ AMPUTEE GUIDES - (10 Pack)        |
| ____ Neoplush socket band                   | ____ Marketing Packets - (10 Pack)     |
| ____ Distal end pads/Per Pkg.               | ____ NEW! MARKETING KIT                |
| ____ Distal end pad/Reticulated SIZE: _____ | ____ PowerPoint Marketing Presentation |



**FLO-TECH® O & P Systems, Inc.**  
 PO Box 462 - 7325 Halseyville Rd  
 Trumansburg, New York 14886

**www.1800FLO-TECH.com**  
**info@1800FLO-TECH.com**

**Phone: (607) 387-3070**  
**Phone: 1-800-FLO-TECH(356-8324)**  
**Fax: (607) 387-3176**

# SYMES Measurement Form

DATE: \_\_\_\_\_ P.O. Number: \_\_\_\_\_

VISA/MASTERCARD ACCT. # \_\_\_\_\_ EXP. DATE: \_\_\_\_\_

NAME AS IT APPEARS ON CARD: \_\_\_\_\_ INVOICE#: \_\_\_\_\_

PATIENT: \_\_\_\_\_ PROSTHETIST: \_\_\_\_\_

COMPANY: \_\_\_\_\_ NEED BY: \_\_\_\_\_

PHONE: \_\_\_\_\_ FAX: \_\_\_\_\_

Shipping Address: \_\_\_\_\_ Billing Address: \_\_\_\_\_

## CUSTOM SYMES ORDER: (PROVIDED BY PRACTITIONER)

LEFT: \_\_\_\_\_ RIGHT: \_\_\_\_\_ From Measurements: \_\_\_\_\_ From Tracing: \_\_\_\_\_ (FAX Two 8" x 11" pieces of paper)

### CIRCLE SOCKETS FOR THIS ORDER:

VCSPS™ OR APOPPS®

#### LENGTHS:

From Mid Patella Tendon To Distal End \_\_\_\_\_ IN \_\_\_\_\_ CM

From Smallest Part of Calf to Distal End \_\_\_\_\_ IN \_\_\_\_\_ CM

From Smallest Part Ankle to Distal End \_\_\_\_\_ IN \_\_\_\_\_ CM

#### AP (ANTERIOR TO POSTERIOR):

At Smallest Part of Ankle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Ball of Symes \_\_\_\_\_ IN \_\_\_\_\_ CM

#### ML (MEDIAL TO LATERAL SIDE):

At Tibial Tubercle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Largest Part of Calf \_\_\_\_\_ IN \_\_\_\_\_ CM

At Smallest Part of Ankle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Ball of Symes \_\_\_\_\_ IN \_\_\_\_\_ CM

#### CIRCUMFERENCES:

At Tibia Tubercle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Widest Part of Calf \_\_\_\_\_ IN \_\_\_\_\_ CM

At Just below Calf Muscle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Mid Point From Distal Calf to Smallest Part of Ankle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Smallest Part of Ankle \_\_\_\_\_ IN \_\_\_\_\_ CM

At Ball of Symes amputation \_\_\_\_\_ IN \_\_\_\_\_ CM

### COMFORT OPTIONS & COMPONENT ITEM PRICES:

| QTY   | DESCRIPTION                                   | QTY   | DESCRIPTION                         |
|-------|---|-------|-------------------------------------|
| _____ | Modification to Standard Socket               | _____ | Each Amputation Sock SIZE: _____    |
| _____ | Distal end pads Per Pkg. (Each) SIZE: _____   | _____ | SS 30 mm Rotatable 4-hole Connector |
| _____ | Waist belt/suspension strap                   | _____ | SS 30 mm SACH Foot Adapter          |
| _____ | Pour & Modify Cast                            | _____ | Foot L _____ R _____ SIZE: _____    |
| _____ | Fork/suspension strap                         | _____ | TI 30 mm Rotatable 4-hole Connector |
| _____ | Velcro® Compatible Neoprene socket band _____ | _____ | TI 30 mm SACH Foot Adapter          |



FLO-TECH® O & P Systems, Inc.  
 PO Box 462 - 7325 Halseyville Rd  
 Trumansburg, New York 14886  
 www.1800FLO-TECH.com

All materials © FLO-TECH® May 11, 1995, Modified August 2011

Phone: 1-800-FLO-TECH (356-8324)  
 Phone: (607) 387-3070  
 Fax: (607) 387-3176  
 E-mail: info@1800FLO-TECH.com